# MEDISANA®

DE **Personenwaage PSA**

GB **Personal scales PSA**

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IT **Pesapersona PSA**

ES **Báscula personal PSA**

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NL **Personenweegschaal PSA**

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SE **Personvåg PSA**

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## IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Should you give this device to another person, it is vital that you also pass on these instructions for use.

### Explanation of symbols

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

#### WARNING

These warning notes must be observed to prevent any injury to the user.

#### CAUTION

These notes must be observed to prevent any damage to the device.

#### NOTE

These notes give you useful additional information on the installation or operation.

Lot number

Manufacturer

## Safety Information

- Any treatment or diet for people who are underweight or overweight requires the qualified advice of a professional (doctor, dietician). The values determined with the scales may be source of useful support.

- ***Risk of tipping over!***

Do not place the scale on an uneven surface.

Never step onto a corner of the scales.

When weighing yourself, do not stand on one side or the edge of the scales.

- ***Risk of slipping!***

Do not stand on the scale with wet feet.

Do not stand on the scale wearing socks.

- Use the device only according to its intended purpose as specified in the instruction manual.
- The warranty will be invalidated if the device is used for purposes other than those for which it is intended.
- The scale has been manufactured for domestic use. It is not suitable for commercial use in hospitals or other medical institutions.
- Do not use the device if it is not working correctly, if it has been dropped or has fallen into water or has been damaged.
- The scales have a measuring range up to 180 kg, 396 lb/28 st. Do not overload the scales.

- Place the scales on a firm, level base. Soft, uneven surfaces are unsuitable for the measuring process and lead to incorrect results.

- Do not use the scale in a location where extremely high temperatures or humidity may occur.

- Keep the scale away from water.

- Handle the scale with care. The scale must not be subjected to impacts or vibration. Do not drop the scale.

- Do not attempt to dismantle the scale, otherwise the warranty will be invalidated. This unit contains no user-maintainable or user-replaceable parts apart from the battery.
- Do not attempt to repair the device yourself in the event of a malfunction since this will invalidate the warranty. Repairs should only be carried out by authorised service centres.
- If the scale is not going to be used for a long period, please remove the batteries.

### Thank you very much

Thank you for your confidence in us and congratulations!

By purchasing the personal scales **PSA**, you have acquired a quality product from **MEDISANA**.

This unit is intended for weighing people.

To ensure that your **MEDISANA** personal scales **PSA** provide successful results for a long time, we recommend that you read carefully through the following notes on their use and care.

### 2.1 Items supplied and packaging

Please check first of all that the unit is complete and is not damaged in any way. If in doubt, do not use the appliance. Send it to a service point.

The following parts are included as standard:

- 1 **MEDISANA PSA** Personal scales
- 3 Type AAA, LR03 1.5 V batteries
- 1 Type AA, LR6 1.5 V battery
- 1 Instructions for use

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

#### **WARNING**

**Please ensure that polybag packing is kept away from reach of children! Risk of suffocation!**

### 2.2 Performance features

- Personal scales with integrated analogue clock
- Large LCD display with luminous white digits on a dark background
- Measuring range up to 180 kg, 396 lbs or 28 st
- 100 g, 0.2 lb graduation
- Switchable between kg, lb and st
- Shatterproof 6 mm safety glass
- Convenient "step-on" switch for weight measurement

### 2.3 Prerequisites for correct measuring results

**Reliable measurements can only be achieved under certain pre-requisites:**

- Always weigh yourself under the same conditions.
- Always weigh yourself at the same time of day.
- Stand up straight and stand still.
- Weigh yourself preferably in the morning, after showering or bathing, and with dry feet.

It is possible to obtain a reliable indication of body weight development from the results if measurements are carried out consistently.

### 3.1 Insert/replace batteries

Before you start up your digital scale, insert the three provided 1.5 V batteries type AAA, LR03 and the 1.5 V battery type AA LR6 (for clock) into the device. To do this, open the battery compartment ❷ beneath the device and insert the batteries. Make sure the polarity is correct (as shown in the battery compartment). Put the lid of the battery compartment back on and push it until it clicks into place.

Change the batteries when the battery change symbol "**Lo**" appears in the display ❶ or when nothing is shown in the display after the device is switched on.

#### **WARNING** **BATTERY SAFETY INFORMATION**

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before putting in the batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of this parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Replace all of the batteries simultaneously!
- Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones!
- Insert the batteries correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge these batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

### 3.2 Change the Unit of Weight

You may have your weight displayed in kilograms (kg), stones (st) or pounds (lb). Use the small switch on the underside of the appliance to set the weight units.

### 3.3 Setting the clock

1. Set the time by turning the small wheel in the battery compartment.
2. If you do not use the clock, remove the battery.

### 3.4 Weighing

1. Place the scale on a stable, firm surface. All four feet must be touching the floor. Uneven and soft flooring affects measuring accuracy.
2. Always weigh yourself on the same scale, at the same place and at the same time. Weigh yourself without clothing, before meals, preferably in the morning before breakfast.
3. If you stand with both feet on the scales, **"0.0 kg"** appears. The values change and the weight unit **"kg"** flashes while your weight is being calculated. When the measurement has stabilised, **"kg"** is shown continuously. Your measured body weight, e.g. **65.1** (65 kg 100 g), will be displayed for approx. 6 seconds. The scales then switch off automatically.

#### NOTE

If the scales have been transported or moved, you can calibrate them by tapping the surface lightly with your foot until **"CAL"** first appears and then **"0.0 kg"** in the display. The scales switch off automatically after approx. 6 seconds.

### 3.5 Error messages

The following error messages may appear in the display:

The scale is overloaded.

The battery is low and must be replaced.

Error during measurement

## 4.1 Care and maintenance

- Remove the batteries before cleaning the unit.
- Never use aggressive cleaning agents or stiff brushes.
- Do not allow any water to get into the unit. Do not immerse the device in water. Do not use the unit again until it is completely dry.
- Do not expose the unit to direct sunlight; protect it against dirt and moisture.
- Do not store the scales upright, as this can cause the batteries to leak.
- Remove the batteries from the unit if it will be unused for a longish period. Otherwise there is a risk that the batteries might leak.

## 4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Please remove the battery before disposing of the equipment. Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.

Consult your municipal authority or your dealer for information about disposal.

## 4.3 Specifications

Name and model	:	<b>MEDISANA PSA</b> Personal scales
Voltage supply	:	scale 4.5 V= , 3 x 1.5 V batteries (type AAA, LR03)
	:	clock 1.5 V= , 1 x 1.5 V battery (type AA, LR6)
Measuring range	:	up to 180 kg, 396 lbs or 28 st
Graduation	:	100 g, 0.2 lb
Autom. switch-off	:	after approx. 6 sec.
Dimensions (L x W x H)	:	approx. 30 x 32 x 4.9 cm
Weight	:	approx. 2.2 kg
Item number	:	40435
EAN number	:	4015588 40435 1

**In accordance with our policy of continual product improvement, we reserve the right for technical and design changes without notice.**

### 5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear (batteries etc.).
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

**MEDISANA AG**

Itterpark 7-9

40724 Hilden

Germany

Tel.: +49 (0) 2103 / 2007-60

Fax: +49 (0) 2103 / 2007-626

eMail: [info@medisana.de](mailto:info@medisana.de)

Internet: [www.medisana.de](http://www.medisana.de)

The service centre address is shown on the last page.

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ALERE Electronics AB